




Transforming Lives®

How to Use the NuStep® T5 Leg Stabilizer

- 1 Insert stability arm into left or right side of seat mount tube on lower front of seat until snap button pops into place.




- 2 Position cuff just back from knee.

Drop strap in between user's legs and wrap strap around leg to secure. Strap does not have to be overly tight. A pre-wrap towel etc. may be used against bare skin.




- 3 Adjust leg alignment by pushing grip forward/up, then slide cuff in or out to set position. There are ten positions (1 – 10) of adjustment.



- 4 Leg Stabilizer should swing approximately 45° during operation (as shown). If it is biased swinging too high or too low outside the 45° range, depress snap button to adjust position up or down. There are three positions (S–M–L) of adjustment.

Verify overall operation by having user start stepping and making any adjustments.





Transforming Lives®

Additional Information

To potentially help move a user's leg into position, install the Leg Stabilizer with the seat swiveled to the side first, and then carefully swivel legs into position.



Two Leg Stabilizers may be used at the same time.

Store Leg Stabilizer(s) on the T5 using a hook on each side. Snap the semi-closed end onto the swivel release lever, and then hang the stabilizer onto the open end.



Or, store the Leg Stabilizer(s) on the optional accessory bag (PN 50704) for neat tidy storage on the seat back.



Additional straps are available for purchase (not shown):

- 51035 Leg Strap, Regular (26"/66cm)
- 51056 Leg Strap, Large (33"/84cm)

Materials:

Nylon, Buna N Rubber, Styrene Butadiene Rubber. Latex Free.

Strap Care Guide:

Machine Wash Cold. No Bleach, Tumble Dry Low.

Product Care Guide:

Use an all-purpose non-abrasive, non-bleach cleaner and a soft cloth.



CAUTION

Leg restraint.

User may not be able to remove leg once installed.

To avoid injury, only align patient's leg to what patient can safely do.

Supervise use of product.